



CUTA Antibody

Product Code	CSB-PA319946ZA01CNE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10951
Immunogen	Recombinant Colletotrichum capsici CUTA protein
Raised In	Rabbit
Species Reactivity	Colletotrichum capsici
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Colletotrichum capsici (Anthracnose fungus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CUTA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.