



aroA Antibody

Product Code	CSB-PA358300ZA01PVU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A2Y4
Immunogen	Recombinant Pseudomonas sp. aroA protein
Raised In	Rabbit
Species Reactivity	Pseudomonas sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pseudomonas sp. (strain PG2982)
Lead Time	Made-to-order (12-14 weeks)
Target Names	aroA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.