



Lingual antimicrobial peptide Antibody

Product Code	CSB-PA389004ZA01BWX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A3RJ36
Immunogen	Recombinant Bubalus bubalis Lingual antimicrobial peptide protein
Raised In	Rabbit
Species Reactivity	Bubalus bubalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Bubalus bubalis (Domestic water buffalo)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Lingual antimicrobial peptide
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.