



fimA Antibody

Product Code	CSB-PA464327ZA01EXZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B2RH54
Immunogen	Recombinant Porphyromonas gingivalis fimA protein
Raised In	Rabbit
Species Reactivity	Porphyromonas gingivalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Porphyromonas gingivalis (strain ATCC 33277 / DSM 20709 / CIP 103683 / JCM 12257 / NCTC 11834 / 2561)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fimA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.