



murl Antibody

Product Code	CSB-PA475624ZA01HUX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B5Z6S2
Immunogen	Recombinant Helicobacter pylori murl protein
Raised In	Rabbit
Species Reactivity	Helicobacter pylori
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Helicobacter pylori (strain G27)
Lead Time	Made-to-order (12-14 weeks)
Target Names	murl
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.