



rpIV Antibody

Product Code	CSB-PA641281ZA01AAAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2GL54
Immunogen	Recombinant Anaplasma phagocytophilum rpIV protein
Raised In	Rabbit
Species Reactivity	Anaplasma phagocytophilum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Anaplasma phagocytophilum (strain HZ)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	rpIV