



Flagellin Antibody

Product Code	CSB-PA23508ZA01GKW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7IZG7
Immunogen	Recombinant Lachnospiraceae bacterium A4 Flagellin protein
Raised In	Rabbit
Species Reactivity	Lachnospiraceae bacterium A4
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Lachnospiraceae bacterium A4
Lead Time	Made-to-order (12-14 weeks)
Target Names	Flagellin
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.