



aro-2 Antibody

Product Code	CSB-PA615522ZA01NHA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12640
Immunogen	Recombinant Neurospora crassa aro-2 protein
Raised In	Rabbit
Species Reactivity	Neurospora crassa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Neurospora crassa (strain ATCC 24698 / 74-OR23-1A / CBS 708.71 / DSM 1257 / FGSC 987)
Lead Time	Made-to-order (12-14 weeks)
Target Names	aro-2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.