



sodB Antibody

Product Code	CSB-PA307060ZA01AJM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83157
Immunogen	Recombinant Anabaena sp. sodB protein
Raised In	Rabbit
Species Reactivity	Anabaena sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Anabaena sp. (strain L31)
Lead Time	Made-to-order (12-14 weeks)
Target Names	sodB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.