



esxA Antibody

Product Code	CSB-PA363657ZA01FSQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A565
Immunogen	Recombinant Mycobacterium bovis esxA protein
Raised In	Rabbit
Species Reactivity	Mycobacterium bovis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
Lead Time	Made-to-order (12-14 weeks)
Target Names	esxA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.