



CLOST_2292 Antibody

Product Code	CSB-PA301302ZA01FUC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80168
Immunogen	Recombinant Acetoanaerobium sticklandii CLOST_2292 protein
Raised In	Rabbit
Species Reactivity	Acetoanaerobium sticklandii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Acetoanaerobium sticklandii (strain ATCC 12662 / DSM 519 / JCM 1433 / NCIMB 10654) (Clostridium sticklandii)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CLOST_2292
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.