



IL17A Antibody

Product Code	CSB-PA011597ZA01RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G1SLF2
Immunogen	Recombinant Oryctolagus cuniculus IL17A protein
Raised In	Rabbit
Species Reactivity	Oryctolagus cuniculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IL17A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.