



# Non-structural glycoprotein 4 Antibody

<b>Product Code</b>	CSB-PA465285ZA01FXW
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P04512
<b>Immunogen</b>	Recombinant Rotavirus A Non-structural glycoprotein 4 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rotavirus A
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugate
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Immunogen Species</b>	Rotavirus A (strain RVA/SA11-Both/G3P5B[2]) (RV-A) (Simian Agent 11 (strain Both))
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	Non-structural glycoprotein 4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.