



glpK Antibody

Product Code	CSB-PA524072ZA01TNR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O66131
Immunogen	Recombinant Thermus thermophilus glpK protein
Raised In	Rabbit
Species Reactivity	Thermus thermophilus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Thermus thermophilus
Lead Time	Made-to-order (12-14 weeks)
Gene Names	glpK