



rpmA Antibody

Product Code	CSB-PA358995ZA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7L8
Immunogen	Recombinant Escherichia coli rpmA protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Escherichia coli (strain K12)
Lead Time	Made-to-order (12-14 weeks)
Target Names	rpmA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.