



CXCL9 Antibody

Product Code	CSB-PA006252ZA01RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	U3KNX2
Immunogen	Recombinant Oryctolagus cuniculus CXCL9 protein
Raised In	Rabbit
Species Reactivity	Oryctolagus cuniculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CXCL9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.