



Polyprotein P1234 Antibody

Product Code	CSB-PA810351ZA01CJAT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8JUX6
Immunogen	Recombinant Chikungunya virus Polyprotein P1234 protein
Raised In	Rabbit
Species Reactivity	Chikungunya virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Chikungunya virus (strain S27-African prototype) (CHIKV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	nsp4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.