



IL11 Antibody

Product Code	CSB-PA324092ZA01MOV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20808
Immunogen	Recombinant Macaca fascicularis IL11 protein
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	IL11