



Poly(3-hydroxybutyrate) depolymerase Antibody

Product Code	CSB-PA319436ZA01FTX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12625
Immunogen	Recombinant <i>Ralstonia pickettii</i> Poly(3-hydroxybutyrate) depolymerase protein
Raised In	Rabbit
Species Reactivity	<i>Ralstonia pickettii</i>
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	<i>Ralstonia pickettii</i> (<i>Burkholderia pickettii</i>)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Poly (3-hydroxybutyrate) depolymerase
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.