



TOP1 Antibody

Product Code	CSB-PA303073ZA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68697
Immunogen	Recombinant Vaccinia virus TOP1 protein
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	TOP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.