



# Recombinant Danio rerio Cellular tumor antigen p53 (tp53)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP024077DIL-B  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P79734   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Danio rerio (Zebrafish) (Brachydanio rerio)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MAQNDSQEFA ELWEKNLIQ PPGGGSCWDI INDEEYLPGS FDPNFFENVL<br>EEQPQPSTLP PTSTVPETSD YPGDHGFRLR FPQSGTAKSV TCTYSPDLNK<br>LFCQLAKTCP VQMVVDVAPP QGSVVRATAI YKKSEHVAEV VRRCPHHERT<br>PDGDNLAPAG HLIRVEGNQR ANYREDNITL RHSVFPYEA PQLGAEWTTV<br>LLNYMCNSSC MGGMNRRLPIL TIITLETQEG QLLGRRSFEV RVCACPGRDR<br>KTEESNFKKD QETKTMAKTT TGTKRSLVKE SSSATLRPEG SKKAKGSSSD<br>EEIFTLQVRG RERYEILKKL NDSLELSDVV PASDAEKYRQ KFMTKNKKEN<br>RESSEPKQGG KLMVKDEGRS DSD |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | tp53   |
| <b>Protein Names</b>     | Recommended name: Cellular tumor antigen p53 Alternative name(s): Tumor suppressor p53   |
| <b>Expression Region</b> | 1-373  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                    |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |