



Replicase large subunit Antibody

Product Code	CSB-PA529603ZA01TOP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O93058
Immunogen	Recombinant Tobacco mosaic virus Replicase large subunit protein
Raised In	Rabbit
Species Reactivity	Tobacco mosaic virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Tobacco mosaic virus (strain OM) (TMV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Replicase large subunit,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.