



hopM1 Antibody

Product Code	CSB-PA803414ZA01FQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q887D0
Immunogen	Recombinant Pseudomonas syringae pv. tomato hopM1 protein
Raised In	Rabbit
Species Reactivity	Pseudomonas syringae pv. tomato
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Pseudomonas syringae pv. tomato (strain ATCC BAA-871 / DC3000)
Lead Time	Made-to-order (12-14 weeks)
Target Names	hopM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.