



Antimicrobial peptide CHP1 Antibody

Product Code	CSB-PA301330ZA01CH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80389
Immunogen	Recombinant Gallus gallus Antimicrobial peptide CHP1 protein
Raised In	Rabbit
Species Reactivity	Gallus gallus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Gallus gallus (Chicken)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Antimicrobial peptide CHP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.