



mxiH Antibody

Product Code	CSB-PA363183ZA01SZB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A223
Immunogen	Recombinant Shigella flexneri mxiH protein
Raised In	Rabbit
Species Reactivity	Shigella flexneri
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Shigella flexneri
Lead Time	Made-to-order (12-14 weeks)
Gene Names	mxiH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.