




# Outer capsid protein VP4 Antibody

<b>Product Code</b>	CSB-PA319263ZA01ROJ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P11196
<b>Immunogen</b>	Recombinant Rotavirus A (strain RVA/Human/United States/DS-1/1976/G2P1B[4]) (RV-A) (Rotavirus A (strain DS1)) Outer capsid protein VP4 protein (247-479aa)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rotavirus A (strain RVA/Human/United States/DS-1/1976/G2P1B[4]) (RV-A) (Rotavirus A (strain DS1))
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Rotavirus A (strain RVA/Human/United States/DS-1/1976/G2P1B[4]) (RV-A) (Rotavirus A (strain DS1))
<b>Target Names</b>	Outer capsid protein VP4 Antibody

## Image



### Western Blot

Positive WB detected in: recombinant protein  
 All lanes: Outer capsid protein VP4 Antibody at 1:1000  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/50000 dilution  
 Predicted band size: 31.4 kDa  
 Observed band size: 32 kDa

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.