



aroA Antibody

Product Code	CSB-PA362185ZA01SXB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07637
Immunogen	Recombinant Salmonella typhimurium aroA protein
Raised In	Rabbit
Species Reactivity	Salmonella typhimurium
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Salmonella typhimurium (strain LT2 / SGSC1412 / ATCC 700720)
Lead Time	Made-to-order (12-14 weeks)
Target Names	aroA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.