



29 kDa outer membrane Antibody

Product Code	CSB-PA306078ZA01ABS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83446
Immunogen	Recombinant Acinetobacter baumannii 29 kDa outer membrane protein
Raised In	Rabbit
Species Reactivity	Acinetobacter baumannii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Acinetobacter baumannii
Lead Time	Made-to-order (12-14 weeks)
Target Names	29 kDa outer membrane protein,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.