



alkB Antibody

Product Code	CSB-PA489400ZA01DQP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B8GWW6
Immunogen	Recombinant Caulobacter vibrioides alkB protein
Raised In	Rabbit
Species Reactivity	Caulobacter vibrioides
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Caulobacter vibrioides (strain NA1000 / CB15N) (Caulobacter crescentus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	alkB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.