



IpxC Antibody

Product Code	CSB-PA460177ZA01AWO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B2I1I9
Immunogen	Recombinant Acinetobacter baumannii IpxC protein
Raised In	Rabbit
Species Reactivity	Acinetobacter baumannii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Acinetobacter baumannii (strain ACICU)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IpxC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.