



Outer capsid protein VP4 Antibody

Product Code	CSB-PA434929ZA01RIR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A9Q1L0
Immunogen	Recombinant Rotavirus X Outer capsid protein VP4 protein
Raised In	Rabbit
Species Reactivity	Rotavirus X
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Rotavirus X (isolate RVX/Human/Bangladesh/NADRV-B219/2002/GXP[X]) (RV ADRV-N) (Rotavirus (isolate novel adult diarrhea rotavirus-B219))
Lead Time	Made-to-order (12-14 weeks)
Target Names	Outer capsid protein VP4,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.