



Dermonecrotic toxin LiSicTox-alphaA1a Antibody

Product Code	CSB-PA317615ZA01LRA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CE80
Immunogen	Recombinant <i>Loxosceles intermedia</i> Dermonecrotic toxin LiSicTox-alphaA1a protein
Raised In	Rabbit
Species Reactivity	<i>Loxosceles intermedia</i>
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	<i>Loxosceles intermedia</i> (Brown spider)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Phospholipase D LiSicTox-alphaA1a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.