



mspA Antibody

Product Code	CSB-PA369785ZA01MVX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0QR29
Immunogen	Recombinant Mycobacterium smegmatis mspA protein
Raised In	Rabbit
Species Reactivity	Mycobacterium smegmatis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Mycobacterium smegmatis (strain ATCC 700084 / mc(2)155)
Lead Time	Made-to-order (12-14 weeks)
Target Names	mspA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.