



pla Antibody

Product Code	CSB-PA325941ZA01YAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17811
Immunogen	Recombinant Yersinia pestis pla protein
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Yersinia pestis
Lead Time	Made-to-order (12-14 weeks)
Gene Names	pla