




KLHDC3 Antibody

Product Code	CSB-PA858414ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQ90
Immunogen	Recombinant Homo sapiens KLHDC3 protein (1-382aa)
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	KLHDC3

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: KLHDC3 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.