



Intermediate capsid protein VP6 Antibody

Product Code	CSB-PA365972ZA01RIX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03530
Immunogen	Recombinant Rotavirus A Intermediate capsid protein VP6 protein
Raised In	Rabbit
Species Reactivity	Rotavirus A
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Rotavirus A (strain RVA/Human/United States/Wa/1974/G1P1A[8]) (RV-A)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Intermediate capsid protein VP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.