



GPC Antibody

Product Code	CSB-PA322925ZA01LCP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17332
Immunogen	Recombinant Lassa virus GPC protein
Raised In	Rabbit
Species Reactivity	Lassa virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Lassa virus (strain GA391) (LASV)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	GPC