



CDA Antibody

Product Code	CSB-PA344641ZA01AJI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50325
Immunogen	Recombinant Amylomyces rouxii CDA protein
Raised In	Rabbit
Species Reactivity	Amylomyces rouxii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Amylomyces rouxii (Filamentous fungus) (Mucor rouxii)
Lead Time	Made-to-order (12-14 weeks)
Target Names	CDA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.