



IpxC Antibody

Product Code	CSB-PA408442ZA01KAX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6T4N9
Immunogen	Recombinant Klebsiella pneumoniae subsp. pneumoniae IpxC protein
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae subsp. pneumoniae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Klebsiella pneumoniae subsp. pneumoniae (strain ATCC 700721 / MGH 78578)
Lead Time	Made-to-order (12-14 weeks)
Target Names	IpxC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.