



ORM2 Antibody

Product Code	CSB-PA017237ZA01GKM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G3QME1
Immunogen	Recombinant Gorilla gorilla gorilla ORM2 protein
Raised In	Rabbit
Species Reactivity	Gorilla gorilla gorilla
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Lead Time	Made-to-order (12-14 weeks)
Target Names	ORM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.