



acpS Antibody

Product Code	CSB-PA431428ZA01ATX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A9NHV3
Immunogen	Recombinant Acholeplasma laidlawii acpS protein
Raised In	Rabbit
Species Reactivity	Acholeplasma laidlawii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Acholeplasma laidlawii (strain PG-8A)
Lead Time	Made-to-order (12-14 weeks)
Target Names	acpS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.