



Collagenase Antibody

Product Code	CSB-PA357515ZA01HRQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08897
Immunogen	Recombinant Hypoderma lineatum Collagenase protein
Raised In	Rabbit
Species Reactivity	Hypoderma lineatum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Hypoderma lineatum (Early cattle grub) (Common cattle grub)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Collagenase
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.