



Silicatein Antibody

Product Code	CSB-PA740986ZA01PGAQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6YD92
Immunogen	Recombinant Petrosia ficiformis Silicatein protein
Raised In	Rabbit
Species Reactivity	Petrosia ficiformis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Petrosia ficiformis (Common Mediterranean sponge)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Silicatein
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.