



PMII Antibody

Product Code	CSB-PA343237ZA01PLO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46925
Immunogen	Recombinant Plasmodium falciparum PMII protein
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Plasmodium falciparum
Lead Time	Made-to-order (12-14 weeks)
Target Names	PMII
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.