



# EPOR Antibody

Product Code	CSB-PA007744ZA01FQF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2KL21
Immunogen	Recombinant Canis lupus familiaris EPOR protein
Raised In	Rabbit
Species Reactivity	Canis lupus familiaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	EPOR