



Genome polyprotein, partial Antibody

Product Code	CSB-PA30568ZA01GKR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A075TBC5
Immunogen	Recombinant Tick-borne encephalitis virus Genome polyprotein,partial protein
Raised In	Rabbit
Species Reactivity	Tick-borne encephalitis virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Tick-borne encephalitis virus
Lead Time	Made-to-order (12-14 weeks)
Target Names	Genome polyprotein,partial
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.