



LOC107913963 Antibody

Product Code	CSB-PA30368ZA01GHB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A1U8K6N0
Immunogen	Recombinant Gossypium hirsutum LOC107913963 protein
Raised In	Rabbit
Species Reactivity	Gossypium hirsutum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Gossypium hirsutum (Upland cotton) (Gossypium mexicanum)
Lead Time	Made-to-order (12-14 weeks)
Target Names	LOC107913963
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.