



Beta-insect excitatory toxin LqqIT1 Antibody

Product Code	CSB-PA324858ZA01LDT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19856
Immunogen	Recombinant Leiurus quinquestriatus quinquestriatus Beta-insect excitatory toxin LqqIT1 protein
Raised In	Rabbit
Species Reactivity	Leiurus quinquestriatus quinquestriatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Leiurus quinquestriatus quinquestriatus (Egyptian scorpion) (Deathstalker scorpion)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	Beta-insect excitatory toxin LqqIT1