



Alpha-like toxin Lqh7 Antibody

Product Code	CSB-PA348581ZA01LDS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59357
Immunogen	Recombinant Leiurus hebraeus Alpha-like toxin Lqh7 protein
Raised In	Rabbit
Species Reactivity	Leiurus hebraeus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Leiurus hebraeus (Deathstalker scorpion) (Leiurus quinquestriatus hebraeus)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Alpha-like toxin Lqh7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.