



cre Antibody

Product Code	CSB-PA362000ZA01FTN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06956
Immunogen	Recombinant Escherichia phage P1 cre protein
Raised In	Rabbit
Species Reactivity	Escherichia phage P1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Immunogen Species	Escherichia phage P1 (Bacteriophage P1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	cre
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.